SOCORE Spinal Osteosynthesis System

Instructions For Use - READ ATTENTIVELY BEFORE USE

The medical device SOCORE is used for the surgical treatment of spinal pathologies. It allows the treatment of spinal deformities, degenerative pathologies, vertebral fractures and tumors. This treatment consists in correcting and stabilizing the spine to help obtaining a fusion between two or more vertebrae and reduce the pain and the disability.

The SOCORE System is intended to assist in the immobilization and stabilization of spinal segments as a complement to fusion in the surgical treatment of spinal pathologies. It helps to reduce pain and disability, with an acceptable rate of rupture and loosening of SOCORE implants.

For best results, a detailed preoperative evaluation, a meticulous surgical technique and adequate postoperative care are mandatory. SOCORE products must be handled and/or implanted by qualified and experienced persons only, who master the surgical procedure, have read the instructions for use and know the device limits. We recommend to the surgeon to do a training course with a surgeon already experimented before starting the use of this technique and this device. It is important that the patient and the surgeon are completely aware of the risks and possible complications associated to this type of surgery.

1) DESCRIPTION AND COMPOSITION OF THE DEVICE

The implants of the SOCORE range are implantable Medical Devices of Class IIb. They are composed of pedicle screws, monoaixial and polyaxial, standard, annulated and reduction; hooks, locking screws, connecting rods, crosslink devices, reinforcement plates and connectors.

The SOCORE Implants are all made of Titanium Alloy according to standard ISO 5832-3. Lumbar short fusion rods are also available in pure Titanium T.40 according to ISO 5832-2, and scoliosis correction rods are available in Cobalt-Chromium alloy according to ISO 5832-12.

It is essential to implant SOCORE products with instrumentation specifically designed for this purpose.

2) INDICATIONS

After meticulous clinical and radiological tests, SOCORE Implants can be used in pediatrics, when the anatomy is favorable, and in adults, by posterior approach, in the following cases:

- | | |
|---|---|
| <ul style="list-style-type: none">- Intervertebral Degenerative Disc Disease- Spinal deformities such as scoliosis, kyphosis and/or lordosis- Spinal trauma | <ul style="list-style-type: none">- Canal stenosis- Spondylolisthesis- Vertebral tumors |
|---|---|

The cervical and occipital system is intended for adults in the cases of spinal trauma and tumors and is used by posterior approach.

3) CONTRA-INDICATIONS

The contra-indications include but are not limited to:

- | | |
|--|---|
| <ul style="list-style-type: none">- Patient not requiring bone graft and fusion- Allergy or intolerance to the materials used- Mental disease, non-cooperative patient- Intense physical activity- Severe osteoporosis- Any case not described in the indications | <ul style="list-style-type: none">- Local or systemic infection- Active smoking- Morbid obesity- Pregnancy- Elderly, severe comorbidity |
|--|---|

4) ADVERSE EFFECTS / POSSIBLE COMPLICATIONS

Potential complications and adverse effects associated with the SOCORE system are identical to those encountered with other spinal systems and **may require a revision surgery**

- | | |
|--|---|
| <ul style="list-style-type: none">- Malposition or mishandling of implants or instruments causing tissue, vascular or nerve damages- Unscrewing, disconnection, torsion and/or migration of any or all implants- Breakage of the implants due to stresses or fatigue- Breakage of bone components (pedicle, vertebra, sacrum)- Skin pressure which may cause irritation, pain, bursitis and potential perforation- Disc degeneration of the adjacent level to the fusion- Neurological injury with loss of function including paralysis (complete or incomplete)- Respiratory disorders (pulmonary embolism, pneumonia, pneumothorax, etc)- Reproductive system compromise, including sterility- Stop of potential growth of the operated segment- Death- The use of these implants in pediatric surgery presents particular risks of bone growth limitation and/or breakage of the implants. In these cases, a close meticulous follow-up is required until the spine is mature. | <ul style="list-style-type: none">- Postoperative change in spinal curvature, loss of correction and/or reduction- Pseudarthrosis- Persistence of pain and pre surgical signs- Infection (superficial or deep)- Dural tears, CSF leakage- Allergy due to implants and materials used- Scattering of wear debris and corrosion- Toxicity due to debris and corrosion, including metallosis, necrosis, tumor, and / or autoimmune disease- Vascular injury, haemorrhage, hematoma- Urinary retention or loss of bladder control- Inability to perform the activities of daily living- Cement extravasation into the neural foramen, spinal canal, soft tissue or intervertebral disc |
|--|---|

5) PREVENTION OF POSTOPERATIVE COMPLICATIONS

- | | |
|---|--|
| <ul style="list-style-type: none">- Avoid extreme forced positions (compression, rotation), avoid all falls- Urgently treat any postoperative infection | <ul style="list-style-type: none">- Avoid all intense physical activities- Avoid excessive body weight- Avoid carrying heavy loads |
|---|--|

6) RISK OF INTERFERENCE WITH MEDICAL IMAGING

MRI / CT: The patient must mention the implants systematically.

The surgeon must also indicate to the patient that in case of investigation and / or treatment (e.g. MRI), he must indicate in advance to the person responsible of the investigation and / or treatment the presence of the medical device in order to limit the reciprocal negative influences.

The SOCORE System has not been evaluated for safety, heating, migration or compatibility in the magnetic resonance environment. MRI imaging quality may be compromised by artifacts if the area of interest is in the exact same area or relatively close to the position of SOCORE devices.

7) PRE-CLEANING / CLEANING

The implants are provided « **NON-STERILE** » and « **for SINGLE USE** »

The implants must be removed from their original packaging, after checking the integrity of this packaging, and checked up to guarantee that they have not been damaged. They must undergo pre-disinfection and cleaning operations with appropriate and compatible products, used in hospitals for this type of implants, before their sterilization. Any cleaning, disinfection and decontamination solution based on aldehyde, alcohol, chlorine, acid or abrasive should be proscribed. The user must validate its pre-disinfection and cleaning procedures according to the applicable standards.

We recommend the exclusive use of mechanized pre-disinfection and cleaning methods with a washer-disinfector (WD).

Automated washing:

Ensure that the WD is in conformity (CE marking), maintained and qualified according to the applicable standards. Use only products compatible with the implants, aldehydes-free, and complying with the Standard Prion Protocol (PSP) if necessary. Follow the instructions for concentration, temperature and duration of action. Use neutral or alkaline pH detergents

Follow the washing procedures to ensure proper protection of the implants.

After cleaning and disinfection in WD, check that all residues have been removed.

If necessary, repeat the cleaning process with brushing until all visible residues are removed.

Recommended Automatic Wash Protocol:

Cycle Step	Minimum time	Recommended temperature	Type of Water / Detergent
Pre-Cleaning	2 minutes	Cold, 20°C	Water
Cleaning	5 minutes	Heated, 55°C	Neutral or Alkaline detergent solution
Neutralization	2 minutes	Cold, 20°C	Water
Rinse	2 minutes	Cold, 20°C	Water
Thermal Rinse	5 minutes	Heated, 93°C	Water
Drying	30 minutes	Air at 99°C	Not applicable (Air)

If necessary, dry residual moisture by means of a complementary drying cycle in the machine or lint-free wiping squares and compressed air.

8) STERILIZATION

Implants are delivered "Non-Sterile" and must be sterilized before use; it is recommended to use a water steam sterilization with a validated cycle including a period of 18 minutes at 134°C and followed by a vacuum drying cycle of 30 minutes minimum. Implants should be prepared so that all surfaces are in direct contact with water vapor. Use absorbent paper between the package and the container and check the absence of moisture inside the container and the package before using the products. Each container and its contents must not exceed 10kg. The implants must be packaged in suitable packaging, handled, stored and transported carefully to guarantee the maintenance of the sterile state until the surgery.

The cycle of sterilization and drying must be validated by the end user by appropriate laboratory techniques and according to the sterilizer manufacturer's recommendations.

Hospitals must qualify and validate their equipment and methods used for steam sterilization, in accordance with current standards for moist heat sterilization of medical devices. The recommendations of the manufacturer of the sterilization equipment must be strictly respected. The sterilization process parameters and the autoclave settings must be monitored regularly.

The hospital assumes final responsibility for the validity of sterilization of products and their maintenance in this state.

Any other method of sterilization (ethylene oxide or oxygen peroxide at low temperature) must not be used and is the sole responsibility of the user.

Sterilization recommendations are given for information purposes only. Under no circumstances can the manufacturer be held responsible for the sterility of sterilized devices within the hospital.

9) REUSE OF DEVICES

The implants that were used by the surgeon, or any other authorized person during the procedure, whether or not installed in the patient, should never be reprocessed or reused. They will be eliminated (see §11.3). Reuse may cause alteration of the claimed performances of the implant as well as contamination of the patient with uncontrolled bioburden.

Implants that have not been used by the surgeon, or any other authorized person during the procedure can be reprocessed, in the absence of signs of wear, damage and loss of performance. The number of reprocessing of unused devices depends on the material of the implants:

For titanium, titanium alloy (ISO 5832-2, ISO 5832-3) and cobalt-chromium alloy (ISO 5832-12) implants, there is no theoretical limit to their reprocessing with the exception of use, damage, wear, or loss of performance. In such cases, reuse involves an alteration of the claimed performances and safety of the implant.

In cases where reprocessing is possible (as described above), the devices must be cleaned and sterilized (see §7 and §8) prior to their use.

10) STORAGE

The implants should be stored carefully in a clean room under normal temperature and humidity conditions.

Implants must be protected from UV rays and all corrosive environments.

11) USE

The surgeon must be experienced with all aspects of the surgical techniques of spinal system in general and the SOCORE system in particular, through his proper surgical training and reading various publications dealing with the spine surgery and the documentation provided. The success of the operation is based on the reduction of pain and disability and the achievement of the bone fusion.

Good patient selection, according to the indications provided, correct follow-up of pre- and postoperative instructions by patients and an adequate surgical procedure are also essential requirements for the success of the operation.

All patients who are considered for implantation of the SOCORE system should be informed about the risks associated to this type of technique. Instructions for the limitation of daily activities and the reduction of stresses on the implants must be provided to the patient. It is an important condition to avoid the failure of the operation and the occurrence of clinical problems.

11.1. Preoperative precautions:

- The surgeon must be completely familiar with all aspects of the device, the instruments and the surgical procedure
- Check before sterilization that an adequate and sufficient range of implants (shapes, sizes) is available
- Check that the implants have not been damaged (absence of scratch, impact, etc...)

- Check before sterilization that all needed instruments for the surgery are available, not worn or distorted and functional
- Check the validity of sterilization

11.2. Intraoperative precautions:

- Proceed according to the surgical procedure provided by the manufacturer
- The correct selection of the type and size of the implants adapted to each other and to the patient as well as their proper use and correct positioning and implantation are extremely important for the success of the operation
- Implant carefully the pedicle screws to avoid incorrect risky positioning
- Correctly insert the locking screw along the axis of implant's head, avoiding tilted position which could damage the implant.
- Check the adequate position and the optimal tightening of all locking screws before closing the soft tissues
- The surgeon must be extremely careful when placing the implants and must pay particular attention to neurological and vascular elements to prevent permanent damage or death
- In revision surgery, all implants have to be replaced with new ones, even if their appearance is correct and do not show signs of wear and fatigue
- Only usable and STERILE implants should be used during the surgical procedure
- During implants manipulation, check their proper surface appearance (absence of scratches, shock, etc.) avoid mechanical stresses during implantation and avoid putting polyaxial screws in extreme position, or any contact or manipulation which can alter their surface or their performance
- Avoid excessive or repetitive bending of the rods and avoid reverse bending
- Do not leave more than one vertebra without implants between two implanted vertebrae
- Components from a different manufacturer or from a different product should not be used in combination with the components of the SOCORE system
- Make bone decortication and use bone graft from the upper to the lower vertebrae of the construct to promote bone fusion
- In case of cement injection, strictly follow the cement manufacturer's instructions for safe and controlled injection. Do not perform bicortical anchorage to prevent leakage of the cement in the retroperitoneal space

11.3. Postoperative precautions:

- Elimination: After pre-disinfection and cleaning, and unless special regulations, the devices to be eliminated will be eliminated according to the methods applied by the hospital, or, if they do not present potential toxic danger, will be sterilized for eventual return to the manufacturer for analysis or expertise.
- The surgeon should give to the patient detailed instructions regarding his activities after surgery, to ensure the survival of the assembly, and should ask him to strictly comply with these instructions (Avoid extreme forced positions, avoid intense physical activity, avoid carrying heavy loads, avoid excessive body weight, avoid stresses on implants, avoid rotation at the spine fusion point, come back for regular postoperative follow-up,...)
- The surgeon will prescribe to the patient an appropriate program of physiotherapy
- It is recommended to make regular postoperative check-ups which allow to check the state of the patient and the good condition of the device
- The respect by the patient of the pre- and postoperative instructions constitutes an essential condition for the success of the operation, especially the reduction of stresses on the implants.
- The implants are intended to stabilize the operative site during the normal healing process and achievement of the fusion. After fusion of the spine, these implants can be removed. The final decision regarding removal of implants depends on the surgeon and the patient. If the device is not removed following completion of its intended use, potential long-term complications may occur, such as corrosion with localized tissue reaction or pain, or unexpected long-term effects such as carcinogenesis.

12) CAUTION


- Implied warranties of merchantability and adequacy to a particular purpose are excluded
- We assume no responsibility in case of combination of our implants with those of another manufacturer
- Implants must not be modified or reworked by the user, except if it complies with SOCORE surgical procedure
- Implantation and eventual removal of implants must be done exclusively with the specific instruments provided by NovaSpine
- Implants cannot support the spine indefinitely without the solid biological support provided by spinal fusion and may develop various defects. These can be a failure at the bone-material interface, rupture of the implants, or bone fracture.
- It has been established that the majority of axial or compression forces are supported by the anterior spinal column. When instrumentation is used to stabilise the spine, adequate anterior support is necessary, either by a surgical procedure (inter somatic reinforcement implants) or by the existing anatomy (satisfactory discal condition). Not maintaining the stability of the anterior column can lead to an overload of the instrumentation and to a failure of implants.
- Despite the care brought by the surgeon to the indications and to the technique of use of the implants, there is a risk of non-fusion, failure, breakage and/or migration of the device, damaging of the bone components and persistence of pain. The surgeon and the patient must be aware of these risks. In such a case, a revision surgery may be necessary, to revise and/or remove devices before a serious injury occurs.

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