SOCORE Sterile Spinal Osteosynthesis System

Instructions For Use (READ ATTENTIVELY BEFORE USE)

The medical device SOCORE is used for the surgical treatment of spinal pathologies. It allows the treatment of spinal deformities, degenerative pathologies, vertebral fractures and tumors. This treatment consists in correcting and stabilizing the spine to promote fusion between two or more vertebrae and reduce the pain and the disability.

The SOCORE System is intended to assist in the immobilization and stabilization of spinal segments as a complement to fusion in the surgical treatment of spinal pathologies. It helps to reduce pain and disability, with an acceptable rate of rupture and loosening of SOCORE implants.

For best results, a detailed preoperative evaluation, a meticulous surgical technique and adequate postoperative care are mandatory. SOCORE products must be handled and/or implanted by qualified and experienced persons only, who master the surgical procedure, have read the instructions for use and know the device limits. We recommend to the surgeon to do a training course with a surgeon already experimented before starting the use of this technique and this device. It is important that the patient and the surgeon are completely aware of the risks and possible complications associated to this type of surgery.

1) DESCRIPTION AND COMPOSITION OF THE DEVICE

The implants of the SOCORE range are implantable Medical Devices of Class IIb. They are composed of pedicle screws, monoaxial and polyaxial, standard, cannulated and reduction; hooks, locking screws, connecting rods, crosslink devices, reinforcement plates and connectors.

The SOCORE Implants are all made of Titanium Alloy according to standard ISO 5832-3. Lumbar short fusion rods are also available in pure Titanium T.40 according to ISO 5832-2.

It is essential to implant SOCORE products with instrumentation specifically designed for this purpose.

2) INDICATIONS

After meticulous clinical and radiological tests, SOCORE Implants can be used in pediatrics and in adults, by posterior and lateral approach, in the following cases:

- | | |
|--|--|
| <ul style="list-style-type: none">- Intervertebral Degenerative Disc Disease- Spinal deformities such as scoliosis, kyphosis and/or lordosis- Spinal trauma- Failure of previous operations | <ul style="list-style-type: none">- Canal stenosis- Spondylolisthesis- Vertebral tumors- Pseudarthrosis |
|--|--|

The cervical and occipital system is intended for adults in the cases of spinal trauma and tumors and is used by posterior approach.

Vertebral staples are used by lateral approach in thoraco-lumbar area.

3) CONTRA-INDICATIONS

The contra-indications include but are not limited to:

- | | |
|--|---|
| <ul style="list-style-type: none">- Patient not requiring bone graft and fusion- Allergy or intolerance to the materials used- Mental disease, non-cooperative patient- Intense physical activity- Any case not described in the indications- Severe osteoporosis | <ul style="list-style-type: none">- Local or systemic infection- Active smoking- Morbid obesity- Pregnancy- Elderly, severe comorbidity |
|--|---|

4) ADVERSE EFFECTS / POSSIBLE COMPLICATIONS

Potential complications and adverse effects associated with the SOCORE system are identical to those encountered with other spinal systems and **may require a revision surgery**

- | | |
|---|--|
| <ul style="list-style-type: none">- Malposition or mishandling of implants or instruments- Damage to soft tissue, vascular or nerve damages- Unscrewing, disconnection, torsion and/or migration of any or all implants- Breakage of the implants due to stresses or fatigue- Breakage of bone components (pedicle, vertebra, sacrum)- Skin pressure which may cause irritation, pain, bursitis and potential perforation- Disc degeneration of the adjacent level to the fusion- Neurological injury with loss of function including paralysis (complete or incomplete)- Respiratory disorders (pulmonary embolism, pneumonia, pneumothorax, etc)- Reproductive system compromise, including sterility- Stop of potential growth of the operated segment- Death | <ul style="list-style-type: none">- Post-operative change in spinal curvature, loss of correction and/or reduction- Pseudarthrosis- Persistence of pain and pre surgical signs- Infection (superficial or deep)- Dural tears, CSF leakage- Allergy due to implants and materials used- Scattering of wear debris and corrosion- Toxicity due to debris and corrosion, including metallosis, necrosis, tumor, and / or autoimmune disease- Vascular injury, haemorrhage, hematoma- Urinary retention or loss of bladder control- Inability to perform the activities of daily living- Cement extravasation into the neural foramen, spinal canal, soft tissue or intervertebral disc |
|---|--|

The use of these implants in pediatric surgery presents particular risks of bone growth limitation and/or breakage of the implants. In these cases, a close meticulous follow-up is required until the spine is mature.

5) PREVENTION OF POSTOPERATIVE COMPLICATIONS

- | | |
|--|--|
| <ul style="list-style-type: none">- Avoid extreme forced positions (compression, rotation), avoid all falls- <u>Urgently treat</u> any post-operative infection | <ul style="list-style-type: none">- Avoid all intense physical activities- Avoid excessive body weight- Avoid carrying heavy loads |
|--|--|

6) RISK OF INTERFERENCE WITH MEDICAL IMAGING

MRI / CT: The patient must mention the implants systematically.

The surgeon must also indicate to the patient that in case of investigation and / or treatment (e.g. MRI), he must indicate in advance to the person responsible of the investigation and / or treatment the presence of the medical device in order to limit the reciprocal negative influences.

The SOCORE System has not been evaluated for safety, warm-up, migration or compatibility in the magnetic resonance environment

7) STERILE PACKAGING

Devices are supplied exclusively in « **STERILE** » form. Packages for each of the components should be intact upon receipt. Once the seal on the sterile package has been broken the product should not be re-sterilized. Damaged packages or products should not be used and should be returned to NovaSpine.

8) RE-USE OF DEVICES

The implants are supplied STERILE and ready to use. The implants that were used by the surgeon, or any other authorized person during the procedure, whether or not installed in the patient, must not be re-used or re-processed. They will be eliminated (see §13). Re-use or re-processing may cause alteration of the claimed performance of the implant as well as contamination of the patient with uncontrolled bioburden.

9) STORAGE

The implants should be stored carefully in a clean room under normal temperature and humidity conditions and protected from UV rays.

Implants must be protected from corrosive environments.

10) USE

The surgeon must be experienced with all aspects of the surgical techniques of spinal system in general and the SOCORE system in particular, through his proper surgical training and reading various publications dealing with the spine surgery and the documentation provided. The success of the operation is based on the reduction of pain and disability and the achievement of the bone fusion.

Good patient selection, according to the indications provided, correct follow-up of pre and post-operative instructions by patients and an adequate surgical procedure are also essential requirements for the success of the operation.

All patients who are considered for implantation of the SOCORE system should be informed about the risks associated to this type of technique. Instructions for the limitation of daily activities and the reduction of stresses on the implants, must be provided to the patient, it is an important condition to avoid the failure of the operation and the occurrence of clinical problems.

11) PRE-OPERATIVE PRECAUTIONS

- The surgeon must be completely familiar with all aspects of the device, the instruments and the surgical procedure
- Check before the surgery that an adequate and sufficient range of implants (shapes, sizes) is available
- Check before sterilization that all needed instruments for the surgery are available, not worn or distorted and functional
- Check the validity of sterilization of the implants
- Check the validity of sterilization of the instruments

12) INTRAOPERATIVE PRECAUTIONS

- Proceed according to the surgical procedure provided by the manufacturer
- The correct selection of type and size of implants adapted to each other and to the patient as well as their proper use and correct implantation and positioning are extremely important for the success of the operation
- Implant carefully the pedicle screws to avoid incorrect risky positioning
- Insert the locking screw in straight position along the axis of implant's head, avoiding tilted position which could damage the implant and tighten firmly on the rod
- Do not use rods with diameter other than Ø5.5 in association with polyaxial pediatric screws (Orbis).
- Check the optimal tightening of all locking screws before closing the soft tissues
- The surgeon must be extremely careful when placing the implants and to pay particular attention to neurological and vascular elements to prevent permanent damage or death
- In revision surgery, all the implants have to be replaced with new ones, even if their appearance is correct and do not show signs of wear and fatigue
- Only STERILE implants should be used during the surgical procedure
- During the manipulation avoid mechanical stresses of the implants during implantation and putting polyaxial screws in extreme position, or any contact or manipulation which can alter their surface or their performance
- Avoid excessive or repetitive bending of the rods and avoid debending at the initial bending point
- Do not leave more than one vertebra without implants between two implanted vertebrae during an assembly
- Components from a different manufacturer or from a different product should not be used in combination with the components of the SOCORE system
- Make bone decortication and use bone graft to achieve a good arthrodesis
- In case of cement injection, strictly follow the cement manufacturer's instructions for safe and controlled injection. Do not take biccortical anchorage to prevent leakage of the cement in the retroperitoneal space

13) POST-OPERATIVE PRECAUTIONS

- Elimination: After pre-disinfection and cleaning, and unless special regulations, the devices to be eliminated will be eliminated according to the methods applied by the hospital, or, if they do not present potential toxic danger, will be sterilized for eventual return to the manufacturer for analysis or expertise.
- The surgeon should give the patient detailed instructions regarding his activities after surgery, to ensure the survival of the assembly, and should ask him to strictly comply with these instructions (Avoid extreme forced positions, avoid intense physical activity, avoid carrying heavy loads, avoid excessive body weight, avoid stresses on implants, avoid turning at the spine fusion point, come back for regular postoperative follow-up,...)
- The surgeon will prescribe to the patient an appropriate program of physiotherapy
- It is recommended to make a regular post-operative check-up which allows to check the state of the patient and the good condition of the device
- The respect by the patient of the pre and post-operative instructions constitutes an essential condition for the success of the operation, especially the reduction of stresses on the implants.

- The implants are intended to stabilize the operative site during the normal healing process and achievement of the fusion. After fusion of the spine, these implants can be removed. The final decision regarding removal of implants depends on the surgeon and the patient. If the device is not removed following completion of its intended use, potential long-term complications may occur, such as corrosion with localized tissue reaction or pain, or unexpected long-term effects such as carcinogenesis.

14) CAUTION

- Implied warranties of merchantability and adequacy to a particular purpose are excluded
- We assume no responsibility if our implants are associated with implants from a different origin
- Implants must not be modified or reworked by the user except if it complies with SOCORE surgical procedure
- Implantation and possible removal of implants must be done exclusively with the specific instruments provided by NovaSpine
- Implants cannot support the spine indefinitely without the solid biological support provided by spinal fusion and may develop various defects, including a failure at the bone-implant interface, rupture or migration of the components, or bone fracture.
- It has been established that the majority of axial or compression forces are supported by the anterior spinal column. When instrumentation is used to stabilise the spine, adequate anterior support is necessary, either by a surgical procedure (inter somatic reinforcement implants) or by the existing anatomy (satisfactory discal condition). The fact of not maintaining the stability of the anterior column can lead to an overload of the instrumentation and to a failing of implants.
- Despite the care brought by the surgeon to the indications and to the technique of use of the implants, there is a risk of non-fusion, failure, breakage and/or migration of the device, of damaging of the bone components and persistence of pain. The surgeon and the patient must be aware of these risks. In such a case, a revision surgery may be necessary, to revise and/or remove devices before a serious injury occurs.

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