DIVA Sterile INTER SOMATIC CAGE

Instructions For Use (READ ATTENTIVELY BEFORE USE)

The DIVA inter somatic cages are used for the surgical treatment of disc pathologies. The treatment consists in distracting, correcting and stabilizing the spine to promote fusion between two or more vertebrae and reduce the pain and the disability.

The DIVA System is intended to help stabilizing spinal segments as a complement to inter-somatic fusion in the surgical treatment of spinal pathologies. It helps to reduce pain and disability, with an acceptable rate of fracture, migration and subsidence of DIVA cages.

For best results, a detailed preoperative evaluation, a meticulous surgical technique and adequate postoperative care are mandatory. DIVA products must be handled and/or implanted by qualified and experienced persons only, who master the surgical procedure, have read the instructions for use and know the device limits. We recommend to the surgeon to do a training course with a surgeon already experimented before starting the use of this technique and this device.

It is important that both the patient and surgeon are fully aware of risks and possible complications associated to this type of surgery.

1) DESCRIPTION AND COMPOSITION OF THE DEVICE

The DIVA inter somatic cages are implantable Medical Devices of Class IIb. They have an anatomical shape and are implanted by impaction into the discal space. One cage is used per segment to fuse except for posterior lumbar approach which requires the use of two cages. Some cervical, anterior and lateral lumbar cages, with a specific design called "Blocked Cages", are fixed to the bone with screws. An additional plate can be added to Non-Blocked cervical, anterior and lateral lumbar cages to prevent their migration. These plates are fixed to the vertebrae with screws. Fusion is made between the two cleaned vertebral end plates using bone grafts previously introduced into the cages, and, when possible, around the cages.

The DIVA inter somatic cages are made of PEEK Optima LT1 in compliance with standard ASTM F2026 and include X-Ray markers made of Titanium alloy according to ISO 5832-3, or Tantalum according to standard ASTM F560. Plates and screws are made of a Titanium alloy, in compliance with ISO 5832-3. It is essential to implant DIVA cages with instrumentation specifically designed for this purpose.

2) INDICATIONS

After meticulous clinical and radiological tests, DIVA cages can be used, for skeletally mature patients, in following cases:

Indications	Lumbar	Cervical
Degenerative disc pathologies (with or without hernia)	✓	✓
Symptomatic loss of discal height	✓	✓
Segmental instability	✓	✓
Cervical spondylotic myelopathy		✓
Low grade degenerative or isthmic spondylolisthesis (with or without reduction)	✓	
Spinal trauma		
Tumoral surgery		
Pseudarthrosis	✓	✓

3) CONTRA-INDICATIONS

The contra-indications include but are not limited to:

- Local or systemic infection
- Allergy or intolerance to the materials used
- Mental disease or non-cooperative patient
- Intense physical activity
- Active smoking
- Patient not requiring bone graft and fusion
- Any case not described in the indications

- Vertebral tumor
- Morbid obesity
- Pregnancy
- Severe osteoporosis
- Vertebral body fracture
- Poor vascular conformation (ALIF)

4) ADVERSE EFFECTS / POSSIBLE COMPLICATIONS

Potential complications and adverse effects associated with the DIVA inter somatic cages system are identical to those encountered with other intersomatic cages systems and **may require a revision surgery**

- Malposition or misusing of the implants or instruments causing tissue, vascular or nerve damages
- Subsidence into vertebral body or retropulsion of cages
- Unscrewing, disassembly, torsion and/or migration of any or all implants
- Breakage of implants due to impaction, stresses or fatigue
- Breakage of bone components (vertebra, sacrum)
- Pressure on nerves, surrounding tissues or organs, ileus
- Disc degeneration of the adjacent level to the fusion
- Neurological injury with loss of function including paralysis (complete or incomplete)
- Vascular injury, haemorrhage, hematoma
- Respiratory disorders (pulmonary embolism, pneumonia, pneumothorax, etc)
- Urinary retention or loss of bladder control
- Inability to perform the activities of daily living
- Death
- Post-operative change in spinal curvature loss of correction and/or reduction
- Dural tears, CSF leakage
- Pseudarthrosis
- Persistence of pain and pre surgical signs
- Infection (superficial or deep)
- Dysphagia and/or perforation of esophagus
- Dysphonia, hoarseness
- Muscle weakness, numbness
- Allergy due to implants and materials used
- Scattering of wear debris and corrosion
- Toxicity due to debris and corrosion, including metallosis, necrosis, tumor and/or autoimmune disease
- Reproductive system compromise, including sterility and retrograde ejaculation (ALIF)

5) PREVENTION OF POSTOPERATIVE COMPLICATIONS

- Avoid extreme forced positions (compression, rotation), avoid all falls
- Avoid all intense physical activities
- Avoid carrying heavy loads
- Avoid excessive body weight
- Urgently treat any post-operative infection

6) RISK OF INTERFERENCE WITH MEDICAL IMAGING

MRI / CT: The patient must mention the implants systematically.

The surgeon must also indicate to the patient that in case of investigation and / or treatment (e.g. MRI), he must indicate in advance to the person responsible of the investigation and / or treatment the presence of the medical device in order to limit the reciprocal negative influences. DIVA System has not been evaluated for been evaluated for safety, warm-up, migration or compatibility in the magnetic resonance environment.

7) STERILE PACKAGING

Devices are supplied exclusively **STERILE**. Packages for each of the components should be intact upon receipt. Once the seal on the sterile package has been broken the product should not be re-sterilized. Damaged packages or products should not be used and should be returned to NovaSpine.

8) RE-USE OF DEVICES

The implants are supplied **STERILE** and ready to use. The implants that were used by the surgeon, or any other authorized person during the procedure, whether or not installed in the patient, must not be re-used or re-processed. They will be eliminated (see §13). Re-use or re-processing may cause alteration of the claimed performance of the implant as well as contamination of the patient with uncontrolled bioburden

9) STORAGE

The implants must be stored with care, in a clean dry room, at normal temperature and protected against UV rays.

Implants must be protected against all corrosive environments.

10) USE

The surgeon must be experienced with all aspects of the surgical techniques of inter somatic cages in general and the DIVA cages in particular, through his proper surgical training and reading various publications dealing with the spine surgery and the documentation provided.

The success of the operation is based on the reduction of pain and disability and the achievement of the bone fusion.

Good patient selection, according to the indications provided, correct follow-up of pre and post-operative instructions by patients and an adequate surgical procedure are also essential requirements for the success of the operation.

All patients who are considered for implantation of the DIVA Cages system should be informed about the risks associated to this type of technique. Instructions for the limitation of daily activities and the reduction of stresses on the implants must be provided to the patient, it is an important condition to avoid the failure of the operation and the occurrence of clinical problems.

11) PRE-OPERATIVE PRECAUTIONS

- The surgeon must be completely familiar with all aspects of the device, the instruments and the surgical procedure
- Check before the surgery that an adequate and sufficient range of implants (shapes, sizes) is available
- Check before sterilization that all needed instruments for the surgery are available, not worn or distorted and functional
- Check the validity of sterilization of the implants
- Check the validity of sterilization of the instruments

12) INTRAOPERATIVE PRECAUTIONS

- Proceed according to the surgical procedure provided by the manufacturer
- The correct selection of the type and size of implants, adapted to each other and to the patient as well as their proper use and correct positioning and implantation are extremely important for the success of the operation.
- Implant carefully the cages, and screws, to avoid incorrect risky positioning
- Avoid strong impaction of Peek cages especially small sizes to avoid their breakage
- The surgeon must be extremely careful when placing the implants and to pay a particular attention to neurological and vascular elements to prevent permanent damage or death.
- In revision surgery, all the implants have to be replaced with new ones, even if their appearance is correct and do not show signs of wear and fatigue
- Only **STERILE** implants should be used during the surgical procedure
- During the manipulation, avoid mechanical stresses of the implants during implantation, or any contact or manipulation which can alter their surface or their performance
- Clean and dissect the vertebral body and use bone graft in the cage and, when possible, around the cages to promote bone fusion

13) POST-OPERATIVE PRECAUTIONS

- Elimination: After pre-disinfection and cleaning, and unless special regulations, the devices to be eliminated will be eliminated according to the methods applied by the hospital, or, if they do not present potential toxic danger, will be sterilized for eventual return to the manufacturer for analysis or expertise.

- The surgeon should give the patient detailed instructions regarding his activities after surgery, to ensure the survival of the assembly, and should ask him to strictly comply with these instructions (Avoid extreme forced positions, avoid intense physical activity, avoid carrying heavy loads, avoid excessive body weight, avoid stresses on the implants, avoid turning at the spine fusion point, come back for regular postoperative follow-up,...)
- The surgeon will prescribe to the patient an appropriate program of physiotherapy
- It is recommended to make a regular post-operative check-up which allows to check the state of the patient and the good resistance of the device
- The respect by the patient of the pre and post-operative instructions constitutes an essential condition for the success of the operation, especially the reduction of stresses on the implants.

14) CAUTION

- Implied warranties of merchantability and adequacy to a particular purpose are excluded
- Implants must not be modified or reworked by the user
- Implantation and possible removal of implants must be done exclusively with the specific instruments provided by NovaSpine
- Implants cannot support the spine indefinitely without the solid biological support provided by spinal fusion and may develop various defects, including a failure at the bone-material interface, rupture or migration of the implants, subsidence in the vertebra or bone fracture.
- **DIVA Intersomatic Non-blocked cages must be used in association with pedicle screws or plates and screws in lumbar use. It is recommended to use plates and screws with non-blocked cervical cages.**
- Despite the care brought by the surgeon to the indications and to the technique of use of the implants, there is a risk of non-fusion, failure, breakage and/or migration of the device, of damaging of the bone components, of subsidence into the vertebral body, and of persistence of pain. The surgeon and the patient must be aware of these risks. In such a case, a revision surgery may be necessary to revise and/or remove devices before a serious injury occurs.

Year of affixing of CE Marking: 2019

CE 1984

**335, Rue Saint Fuscien
80090 Amiens
FRANCE**

**Tel: 33 (0)3 22500731
Fax: 33 (0)3 22463435**

**e-mail:
contact@novaspine.fr**

-  Read the Instructions
-  Consult instructions for use
-  Keep away from sunlight
-  Single Use
-  Do not re-sterilize
-  Keep in Dry Place
-  Do Not Use if Damaged Packaging
-  Sterilized using irradiation
-  Lot number
-  Catalogue number
-  Manufacturer
-  Date of manufacture
-  Use-by date

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