RODD Inter spinous Device

Instructions For Use - READ ATTENTIVELY BEFORE USE

The RODD Interspinous Rotating Device is used for the surgical treatment of lumbar spinal pathologies. The treatment consists in restoring the height between spinous processes and ensuring a dynamic stabilization. It can be used on two adjacent levels maximum.

The RODD system contributes to the reduction of pain and disability, with an acceptable rate of rupture / dislocation and migration of the implant.

For best results, a detailed preoperative evaluation, a meticulous surgical technique and adequate postoperative care are mandatory. RODD products must be handled and/or implanted by qualified and experienced persons only, who master the surgical procedure, have read the instructions of use and know the device limits. Before attempting this technique, surgeons are advised to attend a training course with a surgeon already experienced with the use of the device and this technique. It is important that both the patient and surgeon are fully aware of the risks and possible complications associated to this type of surgery.

1) DESCRIPTION AND COMPOSITION OF THE DEVICE

The RODD Interspinous device is an implantable Medical Device of Class IIb. It is available in different sizes to adapt itself to the spinal anatomy. It is impacted between the spinous processes by posterior approach. One implant is necessary by segment to be treated. Each implant is made of two parts, convex male and concave female, each one is fixed to a spinous process and articulated together. An assembling screw helps avoiding the dislocation of the implant. The two parts of the RODD Interspinous Device are made of PEEK Optima LT1 in compliance with standard ASTM F2026. The assembling screw is made of Titanium alloy according to standard ISO 5832-3.

It is essential to implant the RODD Interspinous device with instrumentation specifically designed for this purpose.

2) INDICATIONS

After meticulous clinical and radiological tests, the RODD Inter spinous device can be used in lumbar area, by posterior approach, for skeletally mature patients, in the following cases:

- Articular facet joint syndrome
- Degenerative or iatrogenic discal deficiency
- Recalibration of degenerative lumbar stenosis
- Symptomatic loss of interspinous height
- Partial arthroctomy
- Recurrence of discal hernia with instability

3) CONTRA-INDICATIONS

The contra-indications include but are not limited to:

- Local or systemic infection
- Vertebral fracture
- Osteoporosis
- Patient needing arthrodesis
- Allergy or intolerance to the materials used
- Mental Disease, Non-cooperative patient
- Active Smoking
- Any case not described in the indications
- Vertebral tumor
- Fever
- Pregnancy
- Morbid obesity
- Severe degenerative stiffness
- Intense physical activity
- Elderly, severe comorbidity

4) ADVERSE EFFECTS / POSSIBLE COMPLICATIONS

Potential complications and adverse effects associated with the RODD Inter Spinous system are identical to those encountered with other inter spinous systems and **may require a revision surgery**:

- Malposition or mishandling of implants or instruments causing tissue, vascular or nerve damages
- Disconnection, torsion and/or migration of any or all implants
- Breakage of the implants due to stresses or fatigue
- Vascular injury, haemorrhage, hematoma
- Neurological injury with loss of function including paralysis (complete or incomplete)
- Respiratory disorders (pulmonary embolism, pneumonia, pneumothorax, etc)
- Inability to perform the activities of daily living
- Death
- Post-operative change in spinal curvature, loss of correction
- Dural tears, CSF leakage
- Persistence of pain and pre surgical signs
- Infection (superficial or deep)
- Breakage of spinous process
- Allergy due to implants and materials used
- Scattering of wear debris and corrosion
- Toxicity due to debris and corrosion, including metallosis, necrosis, tumor and / or autoimmune disease

5) PREVENTION OF POSTOPERATIVE COMPLICATIONS

- Avoid extreme forced positions (compression, extension, rotation), avoid all falls
- Avoid all intense physical activities
- Avoid carrying heavy loads
- Avoid excessive body weight
- Urgently treat any postoperative infection

6) RISK OF INTERFERENCE WITH MEDICAL IMAGING

MRI / CT: The patient must mention the implants systematically.

The surgeon must also indicate to the patient that in case of investigation and / or treatment (e.g. MRI), he must indicate in advance to the person responsible of the investigation and / or treatment the presence of the medical device in order to limit the reciprocal negative influences.

The RODD System has not been evaluated for safety, heating, migration or compatibility in the magnetic resonance environment. MRI imaging quality may be compromised by artifacts if the area of interest is in the exact same area or relatively close to the position of RODD devices.

7) PRE-CLEANING / CLEANING

The implants are provided « **NON-STERILE** » and « **for SINGLE USE** »

The implants must be removed from their original packaging, after checking the integrity of this packaging, and checked up to guarantee that they have not been damaged. They must undergo pre-disinfection and cleaning operations with appropriate and compatible products, used in hospitals for this type of implants, before their sterilization. Any cleaning, disinfection and decontamination solution based on aldehyde, alcohol, chlorine, acid or abrasive should be proscribed. The user must validate its pre-disinfection and cleaning procedures according to the applicable standards.

We recommend the exclusive use of mechanized pre-disinfection and cleaning methods with a washer-disinfector (WD).

7.1. Automated washing:

Ensure that the WD is in conformity (CE marking), maintained and qualified according to the applicable standards.

Use only products compatible with the implants, aldehydes-free, and complying with the Standard Prion Protocol (PSP) if necessary. Follow the instructions for concentration, temperature and duration of action. Use neutral or alkaline pH detergents

Follow the washing procedures to ensure proper protection of the implants.

After cleaning and disinfection in WD, check that all residues have been removed.

If necessary, repeat the cleaning process, with brushing until all visible residues are removed.

Recommended Automatic Wash Protocol:

Cycle Step	Minimum time	Recommended temperature	Type of Water / Detergent
Pre-Cleaning	2 minutes	Cold, 20°C	Water
Cleaning	5 minutes	Heated, 55°C	Neutral or Alkaline detergent solution
Neutralization	2 minutes	Cold, 20°C	Water
Rinse	2 minutes	Cold, 20°C	Water
Thermal Rinse	5 minutes	Heated, 93°C	Water
Drying	30 minutes	Air at 99°C	Not applicable (Air)

If necessary, dry residual moisture by means of a complementary drying cycle in the machine or lint-free wiping squares and compressed air.

8) STERILIZATION

Implants are delivered "Non-Sterile" and **must be sterilized before use**; it is recommended to use a water steam sterilization with a validated cycle including a period of 18 minutes at 134°C and followed by a vacuum drying cycle of 30 minutes minimum. Implants should be prepared so that all surfaces are in direct contact with water vapor. Use absorbent paper between the package and the container and check the absence of moisture inside the container and the package before using the products. Each container and its contents must not exceed 10kg. The implants must be packaged in suitable packaging, handled, stored and transported carefully to guarantee the maintenance of the sterile state until the surgery.

The cycle of sterilization and drying must be validated by the end user by appropriate laboratory techniques and according to the sterilizer manufacturer's recommendations.

Hospitals must qualify and validate their equipment and methods used for steam sterilization, in accordance with current standards for moist heat sterilization of medical devices. The recommendations of the manufacturer of the sterilization equipment must be strictly respected. The sterilization process parameters and the autoclave settings must be monitored regularly.

The hospital assumes final responsibility for the validity of sterilization of products and their maintenance in this state.

Any other method of sterilization (ethylene oxide or oxygen peroxide at low temperature) must not be used and is the sole responsibility of the user.

Sterilization recommendations are given for information purposes only. Under no circumstances can the manufacturer be held responsible for the sterility of sterilized devices within the hospital.

9) REUSE OF DEVICES

The implants that were used by the surgeon, or any other authorized person during the procedure, whether or not installed in the patient, should never be reprocessed or reused. They will be eliminated (see §11.3). Reuse may cause alteration of the claimed performance of the implant as well as contamination of the patient with uncontrolled bioburden.

Implants that have not been used by the surgeon, or any other authorized person during the procedure can be reprocessed, in the absence of signs of wear, damage and loss of performance. The number of reprocessing of unused devices depends on the material of the implants:

For PEEK (ASTM F2026) and titanium alloy (ISO 5832-3) implants, there is no theoretical limit to their reprocessing with the exception of use, damage, wear, or loss of performance. In such cases, reuse involves an alteration of the claimed performance and safety of the implant.

In cases where reprocessing is possible (as described above), the devices must be cleaned and sterilized (see §7 and §8) prior to their reuse.

10) STORAGE

The implants should be stored carefully in a clean room under normal temperature and humidity conditions. Implants must be protected from UV rays and all corrosive environments.

11) USE

The surgeon must be experienced with all aspects of the surgical techniques of interspinous devices in general and the RODD Interspinous device in particular, through his proper surgical training and reading various publications dealing with the spine surgery and the documentation provided. The success of the operation is based on the reduction of pain and disability.

Good patient selection, according to the indications provided, correct follow-up of pre- and postoperative instructions by patients and an adequate surgical procedure are also essential requirements for the success of the operation.

All patients who are considered for implantation of the RODD Interspinous device should be informed about the risks associated to this type of technique. Instructions for the limitation of daily activities and the reduction of stresses on the implants must be provided to the patient. It is an important condition to avoid the failure of the operation and the occurrence of clinical problems.

11.1) Preoperative precautions:

- The surgeon must be completely familiar with all aspects of the device, the instruments and the surgical procedure
- Check before sterilization that an adequate range of implants (shapes, sizes) is available
- Check that the implants have not been damaged (absence of scratch, impact, etc...)
- Check before sterilization that all needed instruments for the surgery are available, not worn or distorted and functional
- Check the validity of sterilization

11.2) Intraoperative precautions:

- Proceed according to the surgical procedure provided by the manufacturer
- The correct selection of the size of the Interspinous device adapted to the patient and the correct handling and positioning of the implant are extremely important for the success of the surgery
- Assemble carefully the two components and screw firmly the assembling screw into the male part
- Impact carefully and as deep as possible the device between the spinous processes to avoid an eventual migration
- The surgeon must be extremely careful when placing the implants and must pay particular attention to neurological and vascular elements to prevent permanent damage
- In revision surgery, all the implants have to be replaced with new ones, even if their appearance is correct and do not show signs of wear and fatigue
- Only usable and STERILE implants should be used during the surgical procedure
- During implants manipulation, check their proper surface appearance (absence of scratches, shock, etc.) avoid mechanical stresses during implantation, or any contact or manipulation which can alter their surface or their performance


11.3) Postoperative precautions:

- Elimination: After pre-disinfection and cleaning, and unless special regulations, the devices to be eliminated will be eliminated according to the methods applied by the hospital, or, if they do not present potential toxic danger, will be sterilized for eventual return to the manufacturer for analysis or expertise.
- The surgeon should give to the patient detailed instructions regarding his activities after surgery, to ensure the survival of the assembly, and should ask him to strictly comply with these instructions (Avoid extreme forced positions, avoid intense physical activity, avoid carrying heavy loads, avoid excessive body weight, avoid stresses on the implants, come back for regular postoperative follow-up,...)
- The surgeon will prescribe to the patient an appropriate program of physiotherapy
- It is recommended to make regular postoperative check-ups which allow to check the state of the patient and the good resistance of the device
- The respect by the patient of the pre- and post-operative instructions constitutes an essential condition for the success of the operation, especially the reduction of stresses on the implants.


12) CAUTION


- Implied warranties of merchantability and adequacy to a particular purpose are excluded
- Implantation and eventual removal of implants must be done with the specific instruments provided by NovaSpine.
- Implants must not be modified or reworked by the user
- Despite the care brought by the surgeon to the indications and to the technique of use of the implants, there is a risk of failure, breakage and/or migration of the device, of breakage of spinous process, and of persistence of pain. The surgeon and the patient must be aware of these risks. In such a case, a revision surgery may be necessary to revise and/or remove devices before a serious injury occurs.

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